

ZIPwire™ Nitinol Hydrophilic Guidewire

Order Number	Diameter	Length	Tip	Taper
Standard Shaft				
M0066802011	0.025 in	150 cm	Straight	3 cm
M0066802051	0.035 in	150 cm	Straight	3 cm
M0066802061	0.035 in	150 cm	Angled	3 cm
M0066802081	0.038 in	150 cm	Straight	3 cm
M0066802091	0.038 in	150 cm	Angled	3 cm
M0066802101	0.038 in	260 cm	Straight	3 cm
Stiff Shaft				
M0066802161	0.025 in	150 cm	Straight	3 cm
M0066802221	0.035 in	150 cm	Straight	3 cm
M0066802231	0.035 in	150 cm	Angled	3 cm
M0066802251	0.038 in	150 cm	Straight	3 cm
Bentson – Type, Standard Shaft				
M0066802141	0.035 in	150 cm	Straight	Bentson 8 cm
M0066802031	0.035 in	150 cm	Angled	Bentson 8 cm
M0066802151	0.038 in	150 cm	Straight	Bentson 8 cm

Quantity: Box 5

Please consult your sales representative for more information and ordering details. To learn more visit www.bostonscientific.eu

1. Torricelli FC, De S, Sarkissian C, et al. Hydrophilic guidewires: evaluation and comparison of their properties and safety. *Urology*. 2013 Nov;82(5):1182-6.

2. Al-Kandari, A. et al. (Eds.) Difficulties in Instrumentation of Endourologic Procedures. Difficult Cases in Endourology. Springer Science & Business Media, 2012.

3. Data on file with Boston Scientific.

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Materials not intended for use in France.

URO-135002-AE Printed in Germany by medicalvision.

ZIPwire™ Nitinol Hydrophilic Guidewire

A smooth hydrophilic guidewire preferred for tortuous anatomy.¹




ZIPwire™

Nitinol Hydrophilic Guidewire

Impacted stones. Tortuous anatomies. Ureteral obstructions.^{1,2}

Medical conditions like these call for a hydrophilic guidewire that's up to the task.¹ Our ZIPwire Nitinol Hydrophilic wire has a kink-resistant core with a proprietary polymer jacket and a lubricious coating to support smooth ureteral access. Its excellent rotational torque³ response enables a smooth response for maximum control. With the flexible hydrophilic design of the ZIPwire Guidewire, you have the power to navigate past obstructions and facilitate successful outcomes – even when faced with challenging cases.^{1,2}



Multiple Tip Options
Available in multiple tip configurations, including the Bentson-type tip designed to provide the additional tip flexibility for negotiation of narrowed ureters and access beyond obstructions

Stiff and Standard Shafts
The Stiff ZIPwire Guidewire is designed to increase pushability, help straighten tortuous anatomy and facilitate advancement of devices²

Highly Radiopaque Jacket
Designed to provide enhanced visibility under fluoroscopy and enable accurate positioning

Torque Vise
May be used for added guidewire steerability and manipulation