




Swiss LithoClast® Trilogy Max

Lithotripter

Enhanced ergonomics. Proven legacy of performance.*

Designed to deliver improved ergonomics without compromising performance, the Swiss LithoClast® Trilogy Max Lithotripter features the enhancements physicians have asked for. It offers a smaller handpiece profile, improved interface and software, and a wider, wireless foot pedal as compared to the Swiss LithoClast Trilogy Lithotripter, its predecessor.

Expect the same efficient performance as the first-generation Trilogy, which was shown to **remove hard stones up to 48%** faster**, on average, than the Olympus ShockPulse-SE Lithotripter, according to bench studies.^{1,2} A prospective, randomized controlled trial also demonstrated that Trilogy had significantly fewer perioperative malfunctions, less suction occlusion, and **49%*** fewer secondary procedures** due to residual stone fragments than the ShockPulse-SE.³



Features and intended benefits when compared to the first-generation Trilogy



Ergonomic handpiece is:

- Designed to be lightweight and comfortable to use in prone and supine patient positioning⁴
- 25% lighter, with a smaller profile⁴
- Fitted with a new cord that reduces the torque effect on the handle by 30%⁴
- Designed to reduce user strain and fatigue due to less torque and weight⁴



Software and interface upgrades include:

- Three presets as well as the ability to customize settings for each case
- An optimized cooling system designed to reduce device downtime due to handpiece overheating during a procedure
- Easier setup for picture-in-picture



Wider, wireless foot pedal is:

- Designed to provide more stability to reduce unintentional pedal movement
- Designed to ease foot placement and activation
- Designed to reduce cord-related hazards and clutter
- Paired automatically
- Powered by a standard AA battery with a low battery indicator designed to minimize case disruption

Ordering information

Order number	Description
M006840202K0	Swiss LithoClast Trilogy Max System Kit
M0068402130	Swiss LithoClast Trilogy Max Handpiece Kit
M0068402150	Swiss LithoClast Trilogy Max Wireless Foot Switch Kit

Please consult your sales representative for more information and ordering details.

* First-generation Trilogy

** Boston Scientific calculation: $(23.79-46.04)/46.04$, 23.79 seconds (Trilogy) vs. 46.04 seconds (ShockPulse-SE), $p<0.01$

*** Boston Scientific calculation: $(17.7-34.7)/34.7$, 17.7% (Trilogy) vs. 34.7% (ShockPulse-SE), $p=0.005$

1. Carlos EC, Wollin DA, Winship BB, et al. In vitro comparison of a novel single probe dual-energy lithotripter to current devices. *J Endourol.* 2018;32:534-540.
2. Bader MJ, Eisel M, Strittmatter F, et al. Comparison of stone elimination capacity and drilling speed of endoscopic clearance lithotripsy devices. *World J Urol.* 2021;39:563-569.
3. Large T, Nottingham CU, Brinkman E, et al. Multi-institutional prospective randomized control trial of novel intracorporeal lithotripters: ShockPulse-SE vs. Trilogy Trial. *J Endourol.* 2021;35:1326-1332.
4. Data on file.

Bench test results may not necessarily be indicative of clinical performance.

Swiss LithoClast® Trilogy Max Lithotripter is manufactured by E.M.S. Electro Medical Systems S.A. and distributed by Boston Scientific Corporation.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

All images are the property of Boston Scientific. All trademarks are the property of their respective owners.

**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

© 2024 Boston Scientific Corporation
or its affiliates. All rights reserved.

URO-1752105-AA JAN 2024