



Maximize the value of the LithoVue™ Ureteroscope

The LithoVue Empower™ Retrieval Deployment Device works with our LithoVue™ Ureteroscope. Now you can turn your trusted LithoVue Ureteroscope into a physician-controlled retrieval device.



Integrate seamlessly

Compatible with Boston Scientific's Zero Tip™, Dakota™, Escape™ and Graspit™ Nitinol Retrieval Devices.

Ordering information

LithoVue Empower™ Retrieval Deployment Device

Product Number	Description	Quantity
M0067919900	LithoVue Empower Retrieval Deployment Device	1 unit

LithoVue™ Single-Use Digital Flexible Ureteroscope

Product Number	Description	Quantity
M0067913500	LithoVue – Global Standard	1 unit
M0067913600	LithoVue – Global Reverse	1 unit
M0067911000*	LithoVue System Workstation Global	1 unit
M0067917000	LithoVue Touch PC	1 unit
M0067914200	LithoVue System Workstation Chart	1 unit
M0067912100	LithoVue System Power Cord – FGGNPSS (France, Greece, Germany, Netherlands, Portugal, Sweden, Spain)	1 unit
M0067912200	LithoVue System Power Cord – UK	1 unit
M0067912300	LithoVue System Power Cord – Italy	1 unit
M0067912400	LithoVue System Power Cord – Denmark	1 unit
M0067912040	LithoVue System Power Cord – Israel	1 unit
M0067912050	LithoVue System Power Cord – Switzerland	1 unit
M0067913480	LithoVue System Replacement Power Adaptor	1 unit
M0067913000	LithoVue System Workstation User's Manual	1 unit
M0067911440	DVI Video Cable	1 unit

* Set Includes: Touch PC, Cart, Power Cord and User's Manual

Want to learn more about the LithoVue Empower Device and how you can take control of stone removal?

Visit www.bostonscientific.eu



Ordering information

Dakota™ Nitinol Stone Retrieval Device with OpenSure™ Handle

Product Number	Description	Quantity
M0063905000	Dakota 1.9F x 8 mm x 120 cm	1 unit
M0063905010	Dakota 1.9F x 11 mm x 120 cm	1 unit

Escape™ Nitinol Stone Retrieval Basket

Product Number	Description	Sheath O.D. (Fr)	Sheath O.D. (mm)	Working Length (cm)	Basket O.D. (mm)
M0063902000	Escape Nitinol Stone Retrieval Basket	1.9	0.63	90	11–15
M0063902010	Escape Nitinol Stone Retrieval Basket	1.9	0.63	120	11–15

Zero Tip™ Nitinol Stone Retrieval Device Basket

Product Number	Description	Sheath O.D. (Fr)	Sheath O.D. (mm)	Working Length (cm)	Basket O.D. (mm)	Sheath Material
M0063901040	Nitinol Zero Tip 4 Wire Basket	1.9	0.63	90	12	Polyimide/PTFE
M0063901050	Nitinol Zero Tip 4 Wire Basket	1.9	0.63	120	12	Polyimide/PTFE
M0063901000	Nitinol Zero Tip 4 Wire Basket	2.4	0.8	90	12	Polyimide/PTFE
M0063901010	Nitinol Zero Tip 4 Wire Basket	2.4	0.8	120	12	Polyimide/PTFE
M0063901020	Nitinol Zero Tip 4 Wire Basket	3.0	1.0	90	16	Polyimide/PTFE
M0063901030	Nitinol Zero Tip 4 Wire Basket	3.0	1.0	120	16	Polyimide/PTFE

Graspit™ Nitinol Stone Retrieval Forceps

Product Number	Description	Sheath O.D. (Fr)	Sheath O.D. (mm)	Working Length (cm)	Basket O.D. (mm)	Sheath Material
M0063204000	Nitinol 2-Prong Forceps	2.6	0.87	90	8	Polyimide/PTFE
M0063204010	Nitinol 2-Prong Forceps	2.6	0.87	120	8	Polyimide/PTFE

1. Joice GA, Ludwig WW, Schwen ZR, et al. Ergonomics and procedure time of novel retrieval deployment device for single surgeon ureteroscopy. Abstract presented at EUS Annual Meeting, May 18, 2018, San Francisco, CA.

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Information not intended for use or distribution in France.

URO-554002-AA JUL 2018 Printed in Germany by medicalvision.



LithoVue Empower™
Retrieval Deployment Device

Take control

Put control in the right hand... or the left.

The LithoVue Empower™ Retrieval Deployment Device takes stone removal procedures to the next level. The critical task of controlling the basket is now in the surgeon's hand while seamlessly integrating with our LithoVue™ Ureteroscope and nitinol baskets.

Discover a revolutionary way to streamline basketing from start to finish.



Keep good form.

The LithoVue Empower Device fits our LithoVue™ Single-Use Digital Flexible Ureteroscope and nitinol baskets without adding muscular workload to your procedure.¹

A better OR experience for everyone.



Communication

Good communication is a cornerstone in any OR. Turning a two-person task into a one-person operation eliminates the chance of miscommunication during stone basketing.



Control

Eliminate the need to coordinate with an assistant during stone retrieval.



Time

Time is a valuable resource – especially in the OR. Now you can retrieve stones single handedly without compromising time.

A bench study demonstrated that stone basketing time with the LithoVue Empower Device was similar to basketing with two surgeons, and was significantly improved as compared to single surgeon basketing.¹

