

Ureteral Catheter Portfolio

BEGIN YOUR PROCEDURES WITH CONFIDENCE

Dual Lumen Ureteral Catheter

By helping deliver irrigation fluids to the ureter and drainage of fluids from the urinary tract, our dual-purpose catheter provides access, advancement or exchange of ancillary devices including guidewires.

- Designed for over-the-wire placement
- Provides two separate working channels for second wire placement or injection

Flexima™ Ureteral Catheter

The single lumen with injection hub provides access and collection.

- Various tip designs to facilitate access to the urinary tract
- Can be placed retrograde
- Radiopaque to aid in visualization

Axxcess™ Ureteral Catheter

The single lumen catheter with injection hub facilitates injection of contrast or establishment of post-procedure drainage collection.

- Tapered tip to provide atraumatic access for delivery of contrast or a guidewire
- Soft Percuflex™ material for ease of placement, graduated markers and radiopaque to facilitate visualization and placement

Imager™ II Ureteral Catheter

The single lumen, torqueable catheter is designed to facilitate difficult access in PCNL and ureteroscopic procedures.

- Biocompatible polymer reinforced with a stainless steel braded wire
- Can be placed retrograde or antegrade



Ordering Information

Dual Lumen Ureteral Catheter

Order Number	UPN	O.D. Size	Working Length	Recommended Guidewire (in)	Injection Lumen	Quantity
405100	M0064051000	10F	54	0.038	0.05	Each

Flexima™ Ureteral Catheter

Order Number	UPN	Description	Catheter O.D. Size	Working Length	Tip O.D. Size	Quantity
400200	M0064002001	Open End	4F	70	-	Each
400201	M0064002011	Open End	5F	70	-	Each
400202	M0064002021	Open End	6F	70	-	Each
400203	M0064002031	Open End	7F	70	-	Each
400204	M0064002041	Open End	8F	70	-	Each
400211	M0064002111	Cone Tip	5F	70	8F	Each
400212	M0064002121	Cone Tip	6F	70	10F	Each
400213	M0064002131	Cone Tip	7F	70	12F	Each
400214	M0064002141	Cone Tip	8F	70	14F	Each
400220	M0064002201	Whistle Tip	4F	70	-	Each
400221	M0064002211	Whistle Tip	5F	70	-	Each
400222	M0064002221	Whistle Tip	6F	70	-	Each
400230	M0064002301	Olive Tip	4F	70	-	Each
400231	M0064002311	Olive Tip	5F	70	-	Each
400236	M0064002361	Spiral Tip	5F	70	-	Each
400237	M0064002371	Spiral Tip	6F	70	-	Each
400241	M0064002411	Spiral Tip	5F	70	14F	Each

Axxcess™ Ureteral Catheter

Order Number	UPN	Description	Catheter O.D. Size	Working Length	Tip O.D. Size	Quantity
400116	M0064001161	Open End	6F	70	-	Box 10

Imager™ II Ureteral Catheter

Order Number	UPN	Description	O.D. Size	Working Length	Quantity
400300	M0064003001	C1	5F	65	Box 5
400301	M0064003011	C2	5F	65	Box 5
400302	M0064003021	Straight	5F	65	Box 5
400303	M0064003031	Bern	5F	65	Box 5
400304	M0064003041	JB1	5F	65	Box 5
400402	M0064004021	Straight	5F	100	Box 5
400403	M0064004031	Bern	5F	100	Box 5
400404	M0064004041	JB1	5F	100	Box 5
400503	M0064005031	Bern	5F	40	Box 5

Please consult your sales representative for more information and ordering details.

To learn more visit www.bostonscientific.com/Urology

Boston Scientific
Advancing science for life™

Caution: U.S. Federal law restricts this device to sale by or on the order of a physician.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. Please check availability with your local sales representative or customer service.

All images are the property of Boston Scientific. All trademarks are the property of their respective owners.

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752
www.bostonscientific.com

© 2019 Boston Scientific Corporation or its affiliates. All rights reserved.

URO-487903-AB SEP 2019

