

Ureteral Catheter Portfolio

BEGIN YOUR PROCEDURES WITH CONFIDENCE

Dual Lumen Ureteral Catheter

By helping deliver irrigation fluids to the ureter and drainage of fluids from the urinary tract, our dual-purpose catheter provides access, advancement or exchange of ancillary devices including guidewires.

- Designed for over-the-wire placement
- Provides two separate working channels for second wire placement or injection

Flexima™ Ureteral Catheter

The single lumen with injection hub provides access and collection.

- Various tip designs to facilitate access to the urinary tract
- Can be placed retrograde
- Radiopaque to aid in visualization

Axxcess™ Ureteral Catheter

The single lumen catheter with injection hub facilitates injection of contrast or establishment of post-procedure drainage collection.

- Tapered tip to provide atraumatic access for delivery of contrast or a guidewire
- Soft Percuflex™ material for ease of placement, graduated markers and radiopaque to facilitate visualization and placement

Imager™ II Ureteral Catheter

The single lumen, torqueable catheter is designed to facilitate access in PCNL and ureteroscopic procedures.

- Biocompatible polymer reinforced with a stainless steel braided wire
- Can be placed retrograde or antegrade



Ordering Information

Dual Lumen Ureteral Catheter

UPN	O.D. Size	Working Length	Recommended Guidewire (in)	Injection Lumen	Quantity
M0064051000	10F	54	0.038	0.05	Each

Flexima™ Ureteral Catheter


UPN	Description	Catheter O.D. Size	Working Length	Tip O.D. Size	Quantity
M0064002001	Open End	4F	70	-	Each
M0064002011	Open End	5F	70	-	Each
M0064002021	Open End	6F	70	-	Each
M0064002031	Open End	7F	70	-	Each
M0064002041	Open End	8F	70	-	Each
M0064002111	Cone Tip	5F	70	8F	Each
M0064002121	Cone Tip	6F	70	10F	Each
M0064002131	Cone Tip	7F	70	12F	Each
M0064002141	Cone Tip	8F	70	14F	Each
M0064002301	Olive Tip	4F	70	-	Each
M0064002311	Olive Tip	5F	70	-	Each

Axxcess™ Ureteral Catheter

UPN	Description	Catheter O.D. Size	Working Length	Tip O.D. Size	Quantity
M0064001161	Open End	6F	70	-	Box 10

Imager™ II Ureteral Catheter

UPN	Description	O.D. Size	Working Length	Quantity
M0064003001	C1	5F	65	Box 5
M0064003011	C2	5F	65	Box 5
M0064003021	Straight	5F	65	Box 5
M0064003031	Bern	5F	65	Box 5
M0064003041	JB1	5F	65	Box 5
M0064004021	Straight	5F	100	Box 5
M0064004031	Bern	5F	100	Box 5
M0064004041	JB1	5F	100	Box 5
M0064005031	Bern	5F	40	Box 5



Please consult your sales representative for more information and ordering details.

To learn more visit www.bostonscientific.eu/Urology

Boston Scientific
Advancing science for life™

www.bostonscientific.eu

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. Please check availability with your local sales representative or customer service.

All images are the property of Boston Scientific. All trademarks are the property of their respective owners.

© 2020 Boston Scientific Corporation or its affiliates. All rights reserved.

URO-487903-AC JUN 2020