

# Ureteral Catheter Portfolio

**BEGIN YOUR PROCEDURES WITH CONFIDENCE**

## Dual Lumen Ureteral Catheter

By helping deliver irrigation fluids to the ureter and drainage of fluids from the urinary tract, our dual-purpose catheter provides access, advancement or exchange of ancillary devices including guidewires.

- Designed for over-the-wire placement
- Provides two separate working channels for second wire placement or injection

## Flexima™ Ureteral Catheter

The single lumen with injection hub provides access and collection.

- Various tip designs to facilitate access to the urinary tract
- Can be placed retrograde
- Radiopaque to aid in visualization

## Axxcess™ Ureteral Catheter

The single lumen catheter with injection hub facilitates injection of contrast or establishment of post-procedure drainage collection.

- Tapered tip to provide atraumatic access for delivery of contrast or a guidewire
- Soft Percuflex™ material for ease of placement, graduated makers and radiopaque to facilitate visualization and placement



## Ordering Information

### Dual Lumen Ureteral Catheter

UPN	O.D. Size	Working Length	Recommended Guidewire (in)	Injection Lumen	Quantity
M0064051000	10F	54	0.038	0.05	Each

### Flexima™ Ureteral Catheter

UPN	Description	Catheter O.D. Size	Working Length	Tip O.D. Size	Quantity
M0064002001	Open End	4F	70	-	Each
M0064002011	Open End	5F	70	-	Each
M0064002021	Open End	6F	70	-	Each
M0064002031	Open End	7F	70	-	Each
M0064002111	Cone Tip	5F	70	8F	Each

### Axxcess™ Ureteral Catheter

UPN	Description	Catheter O.D. Size	Working Length	Tip O.D. Size	Quantity
M0064001161	Open End	6F	70	-	Box 10

Please consult your sales representative for more information and ordering details.

To learn more visit  
[www.bostonscientific.com/EMEA\\_Urology](http://www.bostonscientific.com/EMEA_Urology)

CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. Please check availability with your local sales representative or customer service.

All images are the property of Boston Scientific. All trademarks are the property of their respective owners.

**Boston  
Scientific**  
Advancing science for life™

[www.bostonscientific.eu](http://www.bostonscientific.eu)

© 2018 Boston Scientific Corporation or its affiliates. All rights reserved.  
DINURO2380EA

