


Percuflex™ Plus Ureteral Stent


A firm ureteral stent for maximum control.

Patients with tortuous anatomies and stricture, such as those who are pregnant or undergoing cancer treatment, require ureteral stents firm enough to endure those challenging conditions. Our Percuflex Plus Ureteral Stent with HydroPlus™ Coating offers a strong, smooth surface and a tapered tip that make it easier to navigate through the ureteral orifice. Its high-retention coil helps maintain the pigtail shape and prevent stent migration – making it ideal for long-term stent placement up to 365 days.¹



- **Product name and size to enable quick identification of stent**
- **Graduated circumferential markings every 5cm along the body of the stent**
- **HydroPlus Coating to facilitate advancement up the ureter**
- **Tapered tip design for ease of access through the ureteral orifice**
- **High retention coil strength to help maintain pigtail shape and prevent stent migration**

Material strength and smooth stent surface facilitate navigation of tight strictures




Percuflex Plus Multiple sizes

Available in 4.8F – 8F
with lengths ranging from
10cm – 30cm

Designed for Optimal Drainage


Low profile stent with larger inner lumen for improved ID/OD ratio.

Traditional Stent Percuflex Stent



HydroPlus Hydrophilic Coating

HydroPlus Coating when coated on the Percuflex Stent material becomes lubricious, creating a water-saturated surface for facilitating advancement up the ureter.



Products. Evidence. Solutions.

The advantages keep adding up.

We go beyond delivering proven quality products. Boston Scientific combines a comprehensive portfolio of innovative stents with insightful evidence and intelligent solutions to give you the advantage you need to improve patient care while reducing overall costs.



Percuflex™ Plus Ureteral Stent Ordering Information

Product Code	GTIN Number	Size	Quantity
M0061751950	08714729067566	4.8F x 10cm	Each
M0061751960	08714729067580	4.8F x 12cm	Each
M0061751970	08714729067603	4.8F x 14cm	Each
M0061751980	08714729067627	4.8F x 16cm	Each
M0061751990	08714729067641	4.8F x 18cm	Each
M0061752500	08714729143970	4.8F x 20cm	Each
M0061752510	08714729067719	4.8F x 22cm	Each
M0061752520	08714729056805	4.8F x 24cm	Each
M0061752530	08714729056812	4.8F x 26cm	Each
M0061752540	08714729067757	4.8F x 28cm	Each
M0061752550	08714729067771	4.8F x 30cm	Each
M0061752600	08714729067788	6F x 20cm	Each
M0061752610	08714729067795	6F x 22cm	Each
M0061752620	08714729056829	6F x 24cm	Each
M0061752630	08714729041184	6F x 26cm	Each
M0061752640	08714729041191	6F x 28cm	Each
M0061752650	08714729067849	6F x 30cm	Each
M0061752700	08714729067856	7F x 20cm	Each
M0061752710	08714729067863	7F x 22cm	Each
M0061752720	08714729056836	7F x 24cm	Each
M0061752730	08714729041207	7F x 26cm	Each
M0061752740	08714729041214	7F x 28cm	Each
M0061752750	08714729056843	7F x 30cm	Each
M0061752800	08714729067917	8F x 20cm	Each
M0061752810	08714729067924	8F x 22cm	Each
M0061752820	08714729056850	8F x 24cm	Each
M0061752830	08714729041221	8F x 26cm	Each
M0061752840	08714729067931	8F x 28cm	Each
M0061752850	08714729067948	8F x 30cm	Each

Percuflex Plus Ureteral Stents are also available in various kit configurations with various types of guidewires.

Please consult your sales representative for more information and ordering details.

1. Boston Scientific. Data on file.

All trademarks are the property of their respective owners. Please refer to package insert provided with this product for complete Indications for Use, Contraindications, Warnings, Precautions, Adverse Events, and Instructions prior to using this product. Caution: Federal (U.S.) law restricts this device to sale by or on the order of a physician. Caution: The law restricts this device to sale by or on the order of a physician. Information for the use only in countries with applicable health authority product registrations.

**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752

www.bostonscientific.com

© 2016 Boston Scientific Corporation
or its affiliates. All rights reserved.

URO-373903-AA FEB 2016