


SureDrive™ Steerable Ureteral Stent

The steerable stent
that puts you in control



The faster you can accurately position the right stent in your patient, the better.

The SureDrive™ Steerable Ureteral Stent kit gives you the power to precisely advance, turn and position your Boston Scientific Percuflex™ Plus, Contour™ or Contour VL™ Stent up the ureter and into place.



Advantages of the SureDrive Steerable Stent

- ✓ Easy to advance the stent up the ureter
- ✓ Ability to retract the stent if placed too far up the ureter
- ✓ Ability to rotate the stent for accurate formation of renal coil
- ✓ Ability to position the bladder coil correctly
- ✓ Fast placement of the stent
- ✓ Quick release of the stent



One-handed release system

A secure locking mechanism puts you in complete control during advancement. Once you have achieved the desired placement, a simple quarter turn and push of the thumb knob releases the stent from the positioner – all with a single hand.

With SureDrive™ Steerable Stent, you know with greater certainty that your stent is properly positioned and performing to the highest standards.

Ordering Information

SureDrive Steerable Stent with Percuflex™ Plus

Order Number	Stent Size*
M0061457410	4.8F x 20cm
M0061457420	4.8F x 22cm
M0061457430	4.8F x 24cm
M0061457440	4.8F x 26cm
M0061457450	4.8F x 28cm
M0061457460	6F x 20cm
M0061457470	6F x 22cm
M0061457480	6F x 24cm
M0061457490	6F x 26cm
M0061457500	6F x 28cm
M0061458500	7F x 20cm
M0061458510	7F x 22cm
M0061458520	7F x 24cm
M0061458530	7F x 26cm
M0061458540	7F x 28cm
M0061458550	8F x 22cm
M0061458560	8F x 24cm
M0061458570	8F x 26cm
M0061458580	8F x 28cm

SureDrive Steerable Stent with Contour™

Order Number	Stent Size*
M0061458590	6F x 20cm
M0061458600	6F x 22cm
M0061458610	6F x 24cm
M0061458620	6F x 26cm
M0061458630	6F x 28cm
M0061458640	7F x 20cm
M0061458650	7F x 22cm
M0061458660	7F x 24cm
M0061458670	7F x 26cm
M0061458680	7F x 28cm
M0061458690	8F x 22cm
M0061458700	8F x 24cm
M0061458710	8F x 26cm
M0061458720	8F x 28cm

SureDrive Steerable Stent with Contour VL™

Order Number	Stent Size*
M0061458730	4.8F x 22-30cm
M0061458740	6F x 22-30cm
M0061458750	7F x 22-30cm

* French size x length between pigtails



**Boston
Scientific**
Advancing science for life™

CAUTION: The law restricts this device to sale by or on the order of a physician. Information for the use only in countries with applicable health authority product registrations.

Please refer to package insert provided with this product for complete Indications for Use, Contraindications, Warnings, Precautions, Adverse Events, and Instructions prior to using this product.

All trademarks are the property of their respective owners.

www.bostonscientific.eu

© 2016 by Boston Scientific Corporation
or its affiliates. All rights reserved.

URO-360906-AA JAN 2016