



INTELLANAV STABLEPOINT™

Ablation Catheter



Ordering Information: **INTELLANAV STABLEPOINT** Electrode Configuration: Quadripolar Electrode Spacing: 4-2.5mm - 2.5mm

Model Number	Shaft Size	Tip Size	Curve Style	Shaft Length
M004 RFSDS9620 0	7.5F	7F / 4mm	Standard	110cm
M004 RFSDS9620K2 0	7.5F	7F / 4mm	Large	110cm

Compatible Sheath: 8.5F

SYSTEM COMPONENTS

Model Number	Description
M004 RARC04 0	INTELLANAV STABLEPOINT Catheter Cable
M004 RA6301 US 0	Connection Box - MAESTRO Force Sensing US*
M004 RA6302 0	Connection Box -INTELLAGEN™ STABLEPOINT
M004 RA6501 0	Connection Box -INTELLAGEN™
M004 117 0	METRIQ™ Irrigation Tubing Set*
M004 4100 0	METRIQ Irrigation Pump*
M004 4000 0	MAESTRO 4000™ Controllor*
M004 4010 0	MAESTRO 4000 Pod, 100W*
M004 GEN10020401 0	INTELLAGEN™ RF Generator Kit, which includes the INTELLAGEN™ Foot Pedal
M004 GEN10020402 0	INTELLAGEN™ Irrigation Pump Kit
M004 GEN10020403 0	INTELLAGEN™ Remote Control Kit
M004 RA6501 0	Connection Box - INTELLAGEN™
M004 GEN10020625 1	Tubing Kits

*For customers employing the MAESTRO 4000 Controllor




INTELLANAV STABLEPOINT Catheter Cable




INTELLANAV STABLEPOINT™ Ablation Catheter

Bidirectional Curve Options:

Standard Curve
M004 RFSDS9620 0



Large (K2) Curve
M004 RFSDS9620K2 0



Catheter measurements are approximate.

Catheter configurations are illustrative representations only and may not reflect actual performance.

**Boston
Scientific**
Advancing science for life™

Cardiology
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

Medical Professionals:
1.800.CARDIAC (227.3422)
Customer Service:
1.888.272.1001

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Directions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

©2025 Boston Scientific Corporation
or its affiliates. All rights reserved.
EP-1583306-AB