



THERASPHERE™ Y-90 Glass Microspheres & TRUSELECT™ Microcatheters

Go Further, Treat Better

With the optimal balance of pushability and flexibility, TruSelect Microcatheters can help you map and navigate the most complex vessels and tortuous liver anatomy, allowing you to treat more patients with TheraSphere Y-90 Glass Microspheres.

Designed so you can plan

TruSelect 2.6 flow rates (4.0mL/sec) provide enhanced visibility and flow to support confident mapping, helping you clearly visualize target anatomy and treatment zones.

Engineered to reach further

The 0.021" inner diameter of TruSelect 2.0 enables access to small, distal vessels, allowing you to navigate closer to the tumor.


One platform for every Y-90 case

TruSelect Microcatheters deliver an optimal balance of pushability and flexibility—giving you the right catheter for each step of the procedure without changing platforms or sacrificing trackability.



Actual case image of TruSelect 2.0
delivering TheraSphere

*Image courtesy of
Dr. Travis McKenzie
Intermountain Medical Center
Murray, UT*



THERASPHERE™ Y-90 Glass Microspheres & TRUSELECT™ Microcatheters

TruSelect Microcatheters Ordering Information

TruSelect 2.0 Microcatheter

UPN	Usable Length (cm)	Proximal/distal O.D. (F)	Inner Diameter (in)	Tip Shape	RO Markers	Max Burst Rating (PSI)
M001394101050	105	2.8/2.0	0.021	Straight	1	800
M001394101300	130			Straight	1	
M001394101550	155			Straight	1	
M001394101750	175			Straight	1	
M001394111050	105			Bern	1	
M001394111300	130			Bern	1	
M001394111550	155			Bern	1	
M001394111750	175			Bern	1	

TruSelect 2.6 Microcatheter

UPN	Usable Length (cm)	Proximal/distal O.D. (F)	Inner Diameter (in)	Tip Shape	RO Markers	Max Burst Rating (PSI)
M00139735105010	105	2.9/2.6	0.025	Straight	1	900
M00139735130010	130			Straight	1	
M00139735155010	155			Straight	1	
M00139735155020	155			Straight	2	
M00139735105110	105			Bern	1	
M00139735130110	130			Bern	1	
M00139735155110	155			Bern	1	
M00139735155120	155			Bern	2	



TruSelect™ Microcatheter Indications, Safety, and Warnings

<https://www.bostonscientific.com/truselect-indications>



TruSelect™ 2.6 Microcatheter Indications, Safety, and Warnings

<https://www.bostonscientific.com/truselect-2-6-indications>



TheraSphere™ Y-90 Glass Microspheres Indications, Safety, and Warnings

<https://www.bostonscientific.com/therasphere-indications>

**Boston
Scientific**
Advancing science for life™

Interventional Oncology and Embolization

300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

To order product or for more information
contact customer service at 1.888.272.1001.

© 2026 Boston Scientific Corporation
or its affiliates. All rights reserved.

PI-2460611-AA

Results from case studies are not necessarily predictive of results in other cases. Results in other cases may vary.
All other trademarks are the property of their respective owners.
TheraSphere is a registered trademark of Theragenics Corporation used under license by Boston Scientific Medical
Device Limited, a wholly owned indirect subsidiary of Boston Scientific Corporation.