




TRUSELECT™ 2.6
Microcatheter

Track. Torque. Flow.
All in one.



Next-generation workhorse
that eliminates trade-offs.

Evolve beyond great.



All-in-one performance from a single microcatheter.

The TruSelect™ 2.6 Microcatheter combines the precision of a 2.4 F with the flow rates of a 2.8 F in one tool. Whether you're navigating tortuous anatomy, delivering large embolics, or maintaining high-flow rates, the TruSelect 2.6 Microcatheter offers trackability, torqueability, and flow in a single catheter. With its unique, hand-shapeable tip, it's the right workhorse, every time.



Hand-Shapable Tip

Shapes at room temperature without heat or steam and retains its form once inside the body



Trackability

Smooth, predictable tracking supports confident navigation through tortuous anatomy



Torqueability

Near 1:1 torque response delivers controlled, responsive steering and enhances performance of the hand-shapeable tip



High-Flow Performance

0.025" ID facilitates high-flow performance and compatibility with coils and large embolics

TruSelect 2.6 Microcatheter

UPN	Usable Length (cm)	Proximal/distal O.D. (F)	Inner Diameter (in)	Tip Shape	RO Markers
M00139735105010	105	2.9/2.6	0.025	Straight	1
M00139735130010	130			Straight	1
M00139735155010	155			Straight	1
M00139735155020	155			Straight	2
M00139735105110	105			Bern	1
M00139735130110	130			Bern	1
M00139735155110	155			Bern	1
M00139735155120	155			Bern	2

Boston Scientific
Advancing science for life™

Interventional Oncology and Embolization

300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

To order product or for more information contact customer service at 1.888.272.1001.

© 2026 Boston Scientific Corporation or its affiliates. All rights reserved.

PI-2467307-AA



TruSelect™ 2.6 Microcatheter Indications, Safety, and Warnings

<https://www.bostonscientific.com/truselect-2-6-indications>